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nonacus  
ADVANCING NON-INVASIVE HEALTHCARE

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	Sensitivity	Specificity	NPV
pTa	86%	86%	93%
T1	95%	86%	99%
T2+	89%	86%	97%
G1	76%	86%	96%
G2	92%	86%	97%
G3	92%	86%	95%
NMIBC	89%	86%	92%
MIBC	89%	86%	97%

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
1. Chao-Zhe Zhu, et al, J Cancer. 2019, A review on the accuracy of bladder cancer detection methods.
2. Quentin Alimi, et al, Neurourol Urodyn. 2018. Reliability of urinary cytology and cystoscopy for the screening and diagnosis of bladder cancer in patients with neurogenic bladder: A systematic review.

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3. Ward et al ; Eur Urol Oncol, 2023, Highly Sensitive and Specific Detection of Bladder Cancer via Targeted Ultra-deep Sequencing of Urinary DNA  
4. Ward et al. BJU Int. 2019, Targeted deep sequencing of urothelial bladder cancers and associated urinary DNA: a 23-gene panel with utility for non-invasive diagnosis and risk stratification

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**GALEAS/BLADDER**  
DETECT | REPORT | MONITOR

**GALEAS Bladder**

Patient ID: \_\_\_\_\_ Sample ID: \_\_\_\_\_ Clinician: \_\_\_\_\_ Customer: \_\_\_\_\_  
Patient Name: \_\_\_\_\_ Received: \_\_\_\_\_ Address: \_\_\_\_\_ Address: \_\_\_\_\_  
Date of Birth: \_\_\_\_\_ Processed: \_\_\_\_\_ Report Date: \_\_\_\_\_  
Sample Type: Urine Pellet

**Result Summary: Positive - Variants detected**

Positive test result indicates that, at the time of GALEAS Bladder testing, cancer associated variants were detected and there is a high likelihood that cancer is present.  
Appropriate clinical follow up is required to confirm a clinical diagnosis.

**Variant Details**

Gene	Variant	IMP
TP53	C380T	93.0%
TERT	NA	77.0%

**Test Description**

GALEAS Bladder (this generation on over 770 patient urine samples, has determined a test positive predictive value (PPV) of 87%, negative predictive value (NPV) of 98% and sensitivity of 89% for the detection of all stages of bladder cancer<sup>1,2,3</sup>

The variants in this test have been validated as part of the GALEAS Bladder Triage Hematuria test only. They have not been validated as predictive markers for disease stratification or for informing treatment decisions.

**Positive Explanation:**  
GALEAS Bladder tests for somatic variants in selected regions from across 23 genes. The presence of somatic variants in these regions in urinary DNA has been shown to associate with the presence of bladder cancer. The detection of one or more of somatic variants indicates a high likelihood that cancer is present.

**Negative Explanation:**  
A GALEAS Bladder negative test result, at the time of testing, is determined by the lack of detection of cancer associated somatic variants in the urine sample, suggesting the presence of bladder cancer is unlikely. However, this does not completely exclude the presence of cancer now or in the future.

**QC Status**

PASS

**QC Status Explanation**

There was sufficient read depth across the regions to confidently determine a result.

**Test Limitations**

The test has not been validated as predictive biomarker for disease stratification or for informing treatment decisions.

**Methodology**

DNA was extracted from urine derived cell pellets and collected using the GALEAS Bladder Home Collection Kit. Extracted genomic DNA subsequently underwent target enrichment using the GALEAS Bladder Target Enrichment protocol with sequencing performed using Illumina sequencing by synthesis chemistry.

Data analysis was performed using the GALEAS Bladder analysis pipeline GALEAS Bladder version 23.12.1

Nonacus Limited  
Galeas Bladder (UK) Ltd, 2021-5  
Registered in England  
Birmingham B2 4AG  
United Kingdom  
Register No. 09502779

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